



Rapid 6 Universal Wipes

High-Grade disinfection for hard surfaces



Rapid 6 Universal Wipes are the powerful, water-based solution for the rapid cleaning and disinfection of hard surfaces. Designed as 'extra-wet' wipes to effortlessly provide high grade sporicidal disinfection that doesn't contribute to pathogen resistance in hospitals.

Overview

Using a unique synergistic blend of surfactants, biocides, quaternary ammonium compounds and PHMB, Rapid 6 Universal Wipes provide targeted performance against:

- Gram-positive bacteria
- Gram-negative bacteria
- Fungi
- Mycobacteria
- Enveloped viruses
- Bacterial endospores

Independently tested using 'real-world' conditions and contact times, Rapid 6 Universal Wipes provide 4 hour residual protection against bacteria & spores.

Proven to kill at least 99.9999% of germs in 60 seconds, Rapid 6 Universal Wipes have equal kill rates on both the surface and the wipe itself –

demonstrating that pathogens are being killed and not simply picked up and transferred to adjoining surfaces. This is important in order to minimise the occurrence of bio-burden transfer within hospitals.

Extra-wet, Rapid 6 Universal Wipes help the user achieve even and consistent biocide application - essential in helping reduce the potential for micro-organism cross-resistance to develop.

Rapid 6 Universal Wipes can be used alongside Rapid 6 Foam - simply foam ahead of the wipes to:

- Extend the cleaning surface area of the wipe
- Increase sporicidal performance by up to 100 times

Ready to use and FREE from:

- Alcohol
 - Peracetic Acid
 - Chlorine Dioxide
- Rapid 6 Universal Wipes have no special PPE or storage requirements and are easy for all members of staff to use throughout the hospital.

Test results

Spores

Pathogen	Test	Pass time
<i>Clostridium difficile</i>	EN 17126	4 min.
	EN 13704	1 min.
<i>Bacillus cereus</i>	EN 17126	4 min.
	EN 14347	1 min.
<i>Bacillus subtilis</i>	EN 17126	4 min.
	EN 14347	1 min.

Fungi

<i>Aspergillus niger</i>	EN 1275	5 min.
	EN 1650	1 min.
	EN 13697	15 min.
	EN 13624	5 min.
<i>Candida albicans</i>	EN 14562	15 min.
	EN 1275	1 min.
	EN 1650	5 min.
	EN 13697	1 min.
	EN 13624	1 min.
	EN 14562	15 min.

Mycobacteria

<i>Mycobacterium avium</i>	EN 14348	5 min.
	EN 14563	5 min.
<i>Mycobacterium terrae</i>	EN 14348	5 min.
	EN 14563	5 min.

Bacteria

<i>Enterococcus hirae</i>	EN 1276	1 min.
	EN 13697	5 min.
	EN 13727	1 min.
<i>Escherichia coli</i>	EN 14561	5 min.
	EN 1276	1 min.
	EN 13697	5 min.
<i>Klebsiella pneumoniae</i>	EN 14561	5 min.
	EN 1276	1 min.
	EN 13697	5 min.
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	EN 13727	1 min.
	EN 14561	5 min.
	EN 1040	1 min.
	EN 1276	1 min.
<i>Pseudomonas aeruginosa</i>	EN 13697	5 min.
	EN 13727	1 min.
	EN 14561	5 min.
	EN 1040	1 min.
<i>Staphylococcus aureus</i>	EN 13697	5 min.
	EN 13727	1 min.
	EN 14561	5 min.
	EN 1040	1 min.

Vancomycin-resistant <i>Enterococcus faecalis</i> (VRE)	EN 1276	1 min.
<i>Salmonella enteritidis</i>	EN 1276	1 min.
<i>Salmonella typhimurium</i>	EN 1276	1 min.

Viruses

H1N1 Swine Flu	EN 14476	1 min.
----------------	----------	--------

Other results

4 hour residuality testing against bacteria and spores
Practical "in use" testing against *C. difficile* spores

Contact us for further details on our test results.


RAPID 6 WIPES ARE DESIGNED FOR PROFESSIONAL USE ONLY. USE DISINFECTANTS SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE. ALWAYS FOLLOW CLEANING PROCEDURES AND GUIDELINES.

* For further details, please contact us.

How to use

After performing a risk assessment, wipe from clean to dirty working towards you in straight 'S' shaped lines. Take care not to clean the same area twice. Let area air dry.


- Always use a new wipe for a new surface.
- Don't use the same wipe on multiple surfaces.
- Dispose of wipe if it becomes soiled or dry.



+



Spray Rapid 6 Foam ahead of your wipe to boost cleaning area and performance.



Where to use

Rapid 6 Universal Wipes are designed for use on all hard surfaces (excluding medical devices) and can be used to clean & disinfect amongst others:

- | | | |
|-------------------|---------------------|-------------------------|
| ■ Isolation rooms | ■ Laboratories | ■ Mortuaries |
| ■ Sluice rooms | ■ Blood spills | ■ Bodily fluid spills |
| ■ Macerators | ■ Beds & mattresses | ■ Dressing trolleys |
| ■ Walls & floors | ■ Steel gurneys | ■ Door handles |
| ■ Commodes | ■ Work stations | ■ Toilets, sinks & taps |

Compatibility

Rapid 6 Universal Wipes are compatible with:

- | | |
|----------|-------------|
| ■ Steel | ■ Cotton |
| ■ PVC | ■ Polyester |
| ■ Rubber | ■ Wood |

Pack sizes

Description	Wipe size	Material	Case Size	Order Code
Bucket of 225 wipes	230 x 265mm	31gsm polviscose	Individual	R6-W-225B
Tub of 200 wipes	130 x 130mm	30gsm polviscose	Case of 12	R6-W-200T